



Xpress Micro Dialyzer

MD300


Quick Guide



Quick Guide



1. Preparation. If only one segment is used use, separate it carefully from 8-segmented MD. Don't touch the membrane, please!
2. Buffer Preparation. Pipette dialysis buffer either in a.) a deep well plate $\leq 1,400 \mu\text{l}$, or b.) in a 5 ml-microcentrifuge tube $V \leq 3,400 \mu\text{l}$.
3. Loading the sample. Bring the pipette with sample volume firmly into the round opening. Sample volume should be between 50 and 300 μl .
4. Introduction. Put the MD or the single segment into a) a deep well plate or b) in a microcentrifuge tube as prepared in step 2.



5. Dialysis. One step dialysis can be done in the same microcentrifuge tube or deep well plate. If more than one dialysis step is required, change the position of MD in the deep well plate channels or use a new microcentrifuge tube.
6. Dialysis time. The dialysis time depends on the compound and the cut-off of the semipermeable membrane.
7. Sample retrieval. Set the pipette volume to 330 μl for a sample of 300 μl . Press the pipette button to first stop, hold it, and bring pipette with pipette tip firmly into round opening. Aspirate the sample.
8. Further analysis. Finally, pipette the sample into a microcentrifuge tube or a micro plate.

Manufactured by:


scienovaTM
 innovating for better science

scienova GmbH
 Spitzweidenweg 30
 07743 Jena
 Germany

p: +49 (0) 3641 504 586
 f: +49 (0) 3641 504 587
 e: info@scienova.com
 w: www.scienova.com

© 2022 scienova GmbH | All rights reserved, including pictures and graphics. This print product is protected by copyright. No part of this print product may be reproduced in any form, including pictures and graphics, or utilized by any information storage and retrieval system without written permission from the copyright owner, except for brief quotations embodied in critical articles and reviews.

I 202212