

GeBAflex

Dialysis Tubes



GeBAflex Tubes

Mini (10-250 µl)	6-8kDa 12-14 kDa 25 kDa
Midi (100-800 µl)	1 kDa 3.5 kDa 6-8kDa
Maxi (0.1-3 ml)	3.5 kDa 6-8 kDa 12-14 kDa 25 kDa 50 kDa
Mega (10-20 ml)	1 kDa 3.5 kDa 6-8kDa 12-14 kDa



EASY HANDLING

- Ready-to-use
- Easy handling with standard laboratory pipettes
- Usable in a float



EFFICIENT

- High yield with up to 97% sample recovery
- High-quality membranes made of high-purity regenerated cellulose without heavy metals and sulphides (treated with EDTA)



APPLICATION FIELDS

- Buffer exchange
- Protein Purification
- Sample Concentration
- Formation studies
- Sample preparation of proteins for MALDI-MS