

Mixing Box

Xpress Magnetic Mixing Dialysis Box



Products with Xpress Mixing Box*

- Xpress Micro Dialyzer MD100
- Xpress Micro Dialyzer MD300
- Xpress Equilibrium Dialyzer ED300



EASY HANDLING

- Ready-to-use tool for immediate application on a magnetic stirrer

HANDLING VIDEO



goo.gl/d3wk2C

*product without magnetic stirrer



EFFICIENT

- Higher dialysis efficiency due to buffer circulation during dialysis



APPLICATION FIELDS

- Protein and peptide sample purification
- Protein renaturation
- Re-buffering of protein samples