

# ZelluTrans

## Dialysis Tubes

EN



Type	Flat width (mm)	Diameter (mm)	ml/cm
<b>T1</b> 3.5 kDa	19	12.1	1.2
	46	29.3	6.7
<b>T2</b> 6-8 kDa	25	16	2.0
	32	20.4	3.3
	40	25.5	5.1
	50	31.9	8.0
<b>T3</b> 12-14 kDa	25	15.9	2.0
	33	21	3.5
	45	28.6	6.4
<b>T4</b> 12-14 kDa	25	16	2.0
	45	28.6	6.4



### ADVANTAGES

- Cost-effective dialysis tool with variable volume range
- Very low adsorption of proteins
- Autoclavable



### APPLICATION FIELDS

- Re-buffering of protein samples
- Protein Purification
- Sample Concentration