

ZelluTrans

中文

透析管



类型	平面宽度 (mm)	直径 (mm)	ml/cm
T1 3-5 kDa	19	12.1	1.2
	46	29.3	6.7
	50	31.9	8.0
T2 6-8 kDa	25	16	2.0
	32	20.4	3.3
	40	25.5	5.1
T3 12-14 kDa	25	15.9	2.0
	33	21	3.5
	45	28.6	6.4
T4 12-14 kDa	25	16	2.0
	45	28.6	6.4



优点

- 经济高效的透析工具，容积范围可变
- 蛋白吸附性非常低
- 高温高压消毒



应用范围

- 蛋白样品再缓冲
- 蛋白纯化
- 样品浓缩